

GUIDELINES FOR SEQUENCING SAMPLE SUBMISSION


DNA Type - Genomic DNA; BAC, Plasmid, Cosmid DNA; PCR product

DNA Quantity for library preparation:

- At least 1000 ng for TruSeq library preparation
- Quantitation should be performed using a **PicoGreen** Assay (if you do not have one, we can do it for you)
- ***NanoDrop UV and NanoDrop Fluorometry are NOT compatible with this protocol***
- We highly recommend you to prepare and provide us some extra amount of DNA to ensure timely process of your samples

DNA Quality for library preparation:

- Must be double-stranded
- Should not be a product of whole genome amplification
- Should not be degraded – verified by agarose gel
- ***The DNA sample will NOT be accepted with out the agarose gel image***
- Should contain NO particulate matter
- Should have an OD 260/280 ~ 1.8
- Sample should have a minimal concentration of 50ng/ul in TE, H₂O, or EB
- Sample should NOT have RNA contamination (gDNA or BAC DNA) or primer/primer dimers (PCR product)



Example of gDNA run in 1% agarose gel

RNA Type

The sample RNA Requirements for RNA-seq library preparation:

- total amount of total RNA -1000 ng
- sample volume < 50 μ l
- quantitated by Ribogreen
- pure (OD 260/280 = 2.0)
- DNA free
- quality assessed on an RNA 6000 Nano Chip on the Agilent 2100 Bioanalyzer instrument. A typical RNA sample will produce a smear that ranges from 0.2 kb to 7 kb. This protocol is not designed for preparing small RNA molecules, for example snoRNA, microRNA, tRNA, etc.

For all other types of samples, please contact Alina Akhunova via email akhunova@ksu.edu or by phone (785) 532-1393